

RUBY[®] XL HD

High-Definition Handheld Video Magnifier

Features at a Glance

- 5-inch widescreen LCD
- 2x -14x continuous zoom
- Built-in reading stand
- Easy flip-out handle
- Fast reading with no blur
- Lightweight design at 10.5 oz.
- 3 hours continuous use
- 20 high-contrast color viewing modes
- Freeze Frame with adjustable magnification
- Save up to 80 images and send to computer with USB
- Reading Line and Masks

Why Buy HD?

A high-definition (HD) picture on a handheld magnifier is easier on the eyes at any magnification level, which means less eye fatigue after long periods of use and a better viewing experience, especially when colors and fine details are important.



HD Magnification in an Ultra-Portable Package

The all-new RUBY XL HD puts crystal clear, high-definition magnification in the palm of your hand. At just 10.5 ounces, the slim RUBY XL HD goes anywhere with you. The five-inch screen brings out even the finest details. You can magnify up to 14 times – books, photographs, newspapers, labels, and more. The 20 high-contrast color modes let you adjust text to be easy on your eyes.

Designed for Comfort and Ease of Use

The RUBY XL HD is easy to learn and to use. It snaps open to hold the screen at a natural angle, sliding effortlessly across a document as you read. Color-coded tactile buttons make it easy to adjust magnification, select the color mode, or add a Reading Line or Masks. The convenient fold-out handle has two positions for a firm and balanced grip while reading cans and bottles or when reaching to read items on shelves.

Save Images and Transfer Them to Your PC

A push of a button freezes the image so you can bring it closer for viewing. You can save 80 images – such as phone numbers, serial numbers, and family pictures, and even transfer them to your computer by connecting the RUBY XL HD via USB.



Learn more online at www.FreedomScientific.com

RUBY[®] XL HD

High-Definition Handheld Video Magnifier

Magnification and Camera:

- 2x to 14x continuous zoom
- 5 megapixel autofocus

Standard Color Modes:

- Full color, Black on White, White on Black, Yellow on Blue, Yellow on Black

Configurable Color Modes:

- Black on Yellow, Green on Black, Black on Green, White on Blue, Black on Blue, Blue on Black, Blue on Yellow, Red on Black, Black on Red, Amber on Black, Black on Amber, Violet on Black, Black on Violet, Red on White, Black on Light Blue

Memory:

- 80 images

Battery:

- Lithium-ion polymer rechargeable battery
- 3 hours of continuous use
- 3-hour charge time

Display:

- 5-inch (12.7 cm.) widescreen LCD
- 16:9 aspect ratio

Dimensions:

- 5.44 x 3.45 x 0.78 in.
(13.8 cm. x 8.8 cm. x 2 cm.)

Weight:

- 10.5 oz. (300 g.)

Designed & Assembled
in the USA 



*Read comfortably for hours
using the built-in reading stand
or two-position flip-out handle*



Package Also Includes:

- Carrying case, charger, wrist strap, and USB cable



Call to Learn More
800-444-4443

www.FreedomScientific.com

*Freedom Scientific is the world's leading
manufacturer of assistive technology products for
those who are low vision or blind.*

Visit us online at www.FreedomScientific.com to learn more about these and other innovative products from Freedom Scientific that make living with low vision easier.



PEARL[®]
Portable Reading System for instant
audible access to printed materials



MAGIC[®]
HD Screen Magnification Software
with built-in screen reading features



TOPAZ XL HD[®]
Desktop Video Magnifier with true HD
and the widest available viewing area