


ONYX[®] Portable HD

Portable Video Magnifier



Crisp, Clear Magnification
From Your Desk to the Blackboard

Powerful and Portable HD Magnification!

All of the advanced features you expect -
in a lightweight, compact design!

www.FreedomScientific.com

1-866-569-8474 • 727-803-8000 (outside US)

ONYX[®] Portable HD

Portable Video Magnifier



Displays rich, vibrant colors



Powerful OCR when connected to a PC with OpenBook[®] and GEM[®] software



Folds easily for transport

Features

- ▶ Foldable and easy to carry
- ▶ Easy-to-use remote control
- ▶ Customizable high-contrast colors
- ▶ Save images to an SD Card
- ▶ Connect to a PC or monitor
- ▶ Document, distance, and self-view
- ▶ Scan and read aloud with OpenBook
- ▶ Reading lines, masks, and shades
- ▶ Color-boost contrast

SPECIFICATIONS

Magnification	Weight
1.9x to 141x continuous <i>based on 22-inch monitor</i>	4.3 lbs / 1.95 kg

Dimensions (<i>collapsed</i>)
4.2 in x 13.5 in x 5.3 in / 10.7 cm x 34.3 cm x 13.5 cm

INCLUDED

Carrying case, remote,
USB universal power supply
Video and power cables
GEM[®] software

CONNECTIONS

HDMI
USB 3.0 (video)
USB 2.0 (power)

OPTIONAL

XY reading table

For product details, including advanced options, specifications, and warranty, please visit www.FreedomScientific.com.

750565-001 Rev. B



1-866-569-8474 • 727-803-8000 (outside US)